

User Guide

Version 1.0

1. BEST PRACTISE	2
1.1 Basic Feature: CircRNAs Visualization	2
1.2 Advanced Feature: MRE and RBP sites Visualization on CircRNAs	4
2. HOW TO DOWNLOAD CIRCVIEW AND TEST DATA.....	5
3. HOW TO MANAGE CIRCNRNAS DATA.....	6
3.1 CircRNAs Identification Tool Management.....	6
3.2 CircRNAs Data Management	6
4. HOW TO MANAGE SPECIES DATA.....	7
5. HOW TO FIND CIRCNRNA DATABASES.....	7
6. HOW TO MANAGE MRE AND RBP DATA.....	8
7. HOW TO INSTALL JAVA VIRTUAL MACHINE	9
8. HOW TO INSTALL MYSQL	9
8.1 For Windows.....	9
8.2 For Mac OS.....	9

Last updated 2017/4/25

1. BEST PRACTISE

1.1 Basic Feature: CircRNAs Visualization

1.1.1 Java Virtual Machine needs to be installed before running this program. See **7 HOW TO INSTALL JAVA VIRTUAL MACHINE.**

1.1.2 Download and decompress “CircView.tar.gz” from


<http://github.com/GeneFeng/CircView/blob/master/CircView.tar.gz>

Double click “CircView.jar” to launch the program.


1.1.3 Download and decompress circRNA test data from

<https://github.com/GeneFeng/CircView/blob/master/testdata/human.tar.gz>

1.1.4 Click “CircRNA”->”Load Data” to load circRNA files according species and tools.




1.1.5 Click gene transcript name on list panel to view the image of the circRNAs.



1.1.6 Gene transcript can be searched by name or location.

1.1.7 Click one "Circle" to view details of each circRNA.



1.1.8 Detailed information and image of CircRNAs can be saved for further use.

1.1.9 Click "Analysis"-"Comparison" to make a comparison between circRNAs with different samples and/or tools.

Click "Save as" button to export results to a ".csv" file for further use.

1.2 Advanced Feature: MRE and RBP binding sites Visualization on CircRNAs

1.2.1 MySQL needs to be installed, see **8 HOW TO INSTALL MYSQL**

1.2.2 Restart CircView.jar

1.2.3 Download and decompress MRE data from

http://gb.whu.edu.cn/CircView/testdata/mre_human.tar.gz

1.2.4 Click "MRE"->"Load Data" to load MRE file.

1.2.5 Download and decompress RBP data from

http://gb.whu.edu.cn/CircView/testdata/rbp_human.tar.gz

1.2.6 Click "RBP"->"Load Data" to load RBP file.

1.2.7 Load circRNAs data, see 1.1.4

1.2.8 Check MRE or RBP to add MRE sites (blue lines) or RBP sites (red triangles) to CircRNAs

<http://github.com/GeneFeng/CircView/blob/master/CircView.tar.gz>

Download circRNAs test data from

<https://github.com/GeneFeng/CircView/blob/master/testdata/>

Download MRE data from

http://gb.whu.edu.cn/CircView/testdata/mre_human.tar.gz

Download RBP data from


http://gb.whu.edu.cn/CircView/testdata/rbp_human.tar.gz

3. HOW TO MANAGE CIRCRNAS DATA

3.1 CircRNAs Identification Tool Management

CircView integrates 6 CircRNAs identification tools (CIRCexplorer, circRNA_finder, CIRI, find_circ, Mapsplice, and UROBORUB) by default.


Users can add or delete tool by using menu “CircRNA”->”Add Tool” or “CircRNA”->”Delete Tool”.



3.2 CircRNAs Data Management


CircView can load CircRNAs data directly from output of default 6 CircRNAs identification tools.

Users can also import circRNAs identified by other tools with six tab delimited columns, including chromosome, start position, end position, running number/name, junction reads and strand.




4. HOW TO MANAGE SPECIES DATA

CircView provides 7 species (Human (hg38), Human (hg19), Mouse (mm10), Mouse (mm9), Zebrafish (zv9), Fly (dm6), C.elegans (ce10)) by default.




Users can also add or delete species annotation data with compatible format by using menu “Species”->”Add Species” or “Species”->”Delete Species”.



5. HOW TO FIND CIRCRNA DATABASES

CircView provides links to existing circRNA databases, including Circ2Traits, CircBase, CircInteractome, CircNet, circRNADb, StarBase, and TSCD.




6. HOW TO MANAGE MRE AND RBP DATA

CircRNAs mainly function as sponges for the regulatory elements, such as miRNA response elements (MREs) and RNA binding proteins (RBPs). CircView provides advanced features to display regulatory elements.

This module requires the users to install MySQL locally, see **8 HOW TO INSTALL MYSQL**. Users can load and display the MRE data identified by TargetScan (<http://targetscan.org/>) and/or the RBP data identified by starBase (<http://starbase.sysu.edu.cn/>) or any other software. The format requires five tab delimited columns, including chromosome, start position, end position, MRE/RBP name and description.

Load MRE or RBP file will create table and deposit data into MySQL database, and Clear MRE or RBP will remove data from the database. As the data are persistent, users should not load the same data for more than once.



7. HOW TO INSTALL JAVA VIRTUAL MACHINE

Java Virtual Machine needs to be installed before running this program. Simply access <http://www.java.com>, download Java, and install it.

8. HOW TO INSTALL MYSQL

8.1 For Windows

8.1.1 Download and decompress MySQL Installation file from

http://gb.whu.edu.cn/CircView/MySQL/mysql_windows.tar.gz

8.1.2 Double click “NDP46-KB3045557-x86-x64-AllOS-ENU.exe” to install .NET Framework.

8.1.3 Double click “mysql-installer-community-5.7.16.0.msi” to install MySQL. Please create password “12345” for user root during installation.

8.2 For Mac OS

8.2.1 Download MySQL Installation file from

http://gb.whu.edu.cn/CircView/MySQL/mysql-5.7.17-macos10.12-x86_64.dmg

8.2.2 Double click “mysql-5.7.17-macos10.12-x86_64.dmg” to install MySQL. Please create password “12345” for user root during installation.